## mi**PROTECT**<sup>®</sup> Sampling device A

Ready to use

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
### Sample collection

#### Sample collection from surfaces or solids

(A) Snap the buffer reservoir until the blue little pin inside is broken.

(B) Squeeze the buffer reservoir several times with two fingers to get the buffer to flow into the tube to moisten the swab.

Remove the swab from the tube and skim the swab against the inside of the tube to remove excess liquid from the swab. Do not touch the swab.

(C) Rotate and rub the moistened swab in the sample to collect sample material.

Place the swab with the collected sample into the tube containing the buffer and reassemble both parts.

(D) Rotate the swab several times (10 times) inside the tube to remove the sample material from the swab.

Remove the white cap from the bottom of the sampling device.

#### Sample collection from liquids

Remove the swab from the tube. Do not touch the swab.

(E) Rotate and rub the swab in the liquid sample to collect sample material.

Place the swab with the collected sample into the tube and reassemble both parts.

(A) Snap the buffer reservoir until the blue little pin inside is broken.

(B) Squeeze the buffer reservoir several times with two fingers to get the buffer to flow into the tube to moisten the swab.

(D) Rotate the swab several times (10 times) inside the tube to remove the sample material from the swab.

Remove the white cap from the bottom of the sampling device.